

XTOL Stock - Kodak TMAX 100 Film @ ISO 100 (20C - 68F)

Multitank JOBO - 270mls volume

Chemicals, Prepare -

- Tap Water Rinse - 270mls
- Kodak XTOL Stock - 270mls
- Ilford Ilfostop - 13.5mls + 256.5mls Water
- Ilford Rapid Fix - 54mls + 216mls Water
- Ilford Wetting Agent - 1.5mls + 288.5mls Water
- Distilled Water Rinse - 270mls

JOBO CPP2 - Set Programme to - P & Temp to 20°C

1. **Water Rinse** - 5 minutes.
2. **Dispose**
3. **Developer** - 7 minutes & 15 Seconds.
4. **Dispose**
5. **Stop** - 1 minute.
6. **Dispose**
7. **Fixer** - 5 minutes.
8. **Dispose**
9. **Water Rinse** - 1 minute.
10. **Dispose**
11. **Water Rinse** - 2 minute.
12. **Dispose**
13. **Water Rinse** - 2 minute.
14. **Dispose**
15. **Distilled Water Rinse** - 5 minute.
16. **Dispose**
17. **Transfer to Open Container and remove from Spool.**
18. **Wetting Agent** - Agitate for 1 minute.